



ARABIC ENGINEERING OFFICE  
المكتب الهندسي العربي  
EXCELLENCE SINCE 1977

## Preventive maintenance checklist

Department name:

Device name:

SW Version:

Model:

SN:

| Preparations        |                                                                    |     |    |     |
|---------------------|--------------------------------------------------------------------|-----|----|-----|
|                     | Task                                                               | Yes | No | N/A |
| 1                   | Performed general inspection/visual check of the system.           |     |    |     |
| 2                   | Checked gas supply hoses and connections.                          |     |    |     |
| 3                   | Checked system serial number, software version and operating time. |     |    |     |
| 4                   | Checked technical log and test result log.                         |     |    |     |
| 5                   | Performed pre-use check.                                           |     |    |     |
| Parts to be checked |                                                                    |     |    |     |
|                     | Task                                                               | Yes | No | N/A |
| 6                   | Checked power cable.                                               |     |    |     |
| 7                   | Checked control cable.                                             |     |    |     |
| 8                   | Checked user interface.                                            |     |    |     |
| 9                   | Checked air & oxygen gas supply pressure.                          |     |    |     |
| Parts to be cleaned |                                                                    |     |    |     |
|                     | Task                                                               | Yes | No | N/A |
| 10                  | Expiratory valve membrane.                                         |     |    |     |
| 11                  | Safety valve membrane.                                             |     |    |     |
| 12                  | Fan filter.                                                        |     |    |     |

ملاحظات القائم بالصيانة :

توقيع ممثل المستشفى

تاريخ وتوقيع القائم بالصيانة